

## Plusfactor™ Reagent

Plusfactor™ Cat. No. AB-PF-0001 Size: 1 ml

Plusfactor™ Cat. No. AB-PF-0004 Size: 4 x 1 ml

Store at +4°C (do not freeze)

## Description

Lipofector EXT™ Reagents are packaged with Lipofector™ Reagent (Cat. No. AB-LF-Q001) and Plusfactor™ Reagent (Cat. No. AB-PF-0001).

Lipofector™ Reagent (Cat. No. AB-LF-Q001) is designed for the transfection of DNA or RNA into eukaryotic cells. Lipofector™ is a polycationic liposomal reagent, offering outstanding transfection efficiency.

Transfection activity can be enhanced by using Plusfactor™ Reagent (Cat. No. AB-PF-0001) to pre-complex the DNA or RNA

## Guidelines for Transfection

1. Prepare complexes using the amount of DNA or RNA, Plusfactor™ and Lipofector™ recommended on below table. Optimizations are necessary. We recommend using serum free medium to dilute Lipofector™, Plusfactor™ and DNA or RNA.

To transfect cells in different tissue culture formats, vary the amounts of Lipofector™, Plusfactor™, DNA or RNA, cells, and medium as shown in the table 1 and table 2

**Table 1. Reagent quantities for general recommendation**

Culture vessel	Surface area (cm <sup>2</sup> )	Vol. of plating medium	DNA or RNA (µg) in media vol. (µl)	Plusfactor™ (µl) in media vol. (µl)	Lipofector™ (µl) in media vol. (µl)	Medium exchange vol.
96-well	0.3	100 µl	0.1 µg in 25 µl	0.3 µl in 25 µl	0.25 µl in 25 µl	60 µl
48-well	0.7	200 µl	0.2 µg in 50 µl	0.6 µl in 50 µl	0.5 µl in 50 µl	120 µl
24-well	2	500 µl	0.4 µg in 50 µl	1.2 µl in 50 µl	1 µl in 50 µl	300 µl
12-well	4	1 ml	0.8 µg in 50 µl	2.4 µl in 50 µl	2 µl in 50 µl	600 µl
6-well	10	2.5 ml	1.5 µg in 100 µl	4.5 µl in 100 µl	4 µl in 100 µl	1.5 ml
60-mm	20	10 ml	3 µg in 100 µl	9 µl in 100 µl	8 µl in 100 µl	6 ml
100-mm	56	30 ml	6 µg in 200 µl	18 µl in 200 µl	20 µl in 200 µl	20 ml

**Table 2. . Reagent quantities for optimizing transfections**

<b>Cells</b>	<b>DNA or RNA</b>	<b>Plusfactor™</b>	<b>Lipofector™</b>
Censitive cells <b>Note 1</b>	0.25 ug	0.5 µl – 1.0 µl	0.5 µl – 1.25 µl
Most cell lines	0.5 ug	1 µl – 2.5µl	1 µl – 3.0 µl
	0.75 ug	1.5 µl – 3.75 µl	1.5 µl – 4.5 µl
Suspension and robust cells <b>Note 2</b>	1 ug	2 µl – 5 µl	2 µl – 6 µl

**Note 1.** (Examples are HT1080 and Hela)

**Note 2.** (Examples are MCF7, Jurkat, HL60 and A549)

2. Don't add antibiotics to media during transfection procedure.
3. 70-90% confluence at the time of transfection is recommended for high efficiency and to minimize cytotoxicity. Optimization should be necessary.
4. Test serum-free media for compatibility with Lipofector™ since some serum-free formulations may inhibit liposomal transfection.
5. To avoid microbial contamination all solutions should be sterile-filtered before use and subsequently be handled under aseptic conditions, as is common practice for handling cell cultures.

### **Transfection Procedure**

Use the following procedure to transfect adherent mammalian cells in a 24- well format. For other formats, see table 1

1. **[Cell culture]** Plate the cells the day before the transfection experiment. The appropriate plating density for a particular cell line will depend on the growth rate and the shape of the cells. The cells should be 70-90% confluent on the day of transfection. As a general guideline, plate  $2-6 \times 10^4$  cells in 500 µl culture medium with the usual amount of serum.
2. **[Medium exchange]** While complexes are forming, remove the growth medium from cells and replace with 250 ul transfection medium without serum (usually the cell growth medium without serum.).
3. **[DNA or RNA dilution]** Dilute 0.4 µg (or 0.25 – 1.0 ug) DNA or RNA in 50 µl of serum free medium (or other appropriate medium) without serum and mix gently.

4. **[Plusfactor dilution]** Dilute 1.2  $\mu$ l (or 0.5 – 5  $\mu$ l) Plusfactor™ in 50  $\mu$ l of serum free medium (or other appropriate medium) without serum and mix gently.

5. **[Premixing]** Combine the diluted DNA (or RNA) (from step3) and diluted Plusfactor™ (from step 4). Mix gently and incubate for 10-15 minutes at room Temperature or 4 °C. (approximate total volume = 100  $\mu$ l)

**Note 3.** Some plasmid or high density of DNA or RNA are better in 4 °C incubation

6. **[Lipofactor dilution]** Dilute 1  $\mu$ l (or 0.5 – 6  $\mu$ l) Lipofactor™ in 50  $\mu$ l of AB-OTM™ optimized transfection medium (or other appropriate medium) without serum and mix gently.

7. **[Complexes formation]** Combine the diluted DNA (or RNA) / Plusfactor™ (from step 5) and diluted Lipofactor™ (from step 6). Mix gently and incubate for 15 - 30 minutes at room Temperature or 4 °C. (approximate total volume = 150  $\mu$ l)

8. **[Transfection]** Add the 150  $\mu$ l of diluted complexes (DNA or RNA + Plusfactor™ + Lipofactor™ / from step 7) to each well. Mix gently by rocking the plate.

9. **[Cell culture]** Incubate cells at 37°C in a CO<sub>2</sub> incubator for 4 - 6 hours. Add 400  $\mu$ l of growth medium containing 2X the normal concentration of serum without removing the transfection mixture. If toxicity is problem, replace medium with fresh, complete medium (with normal amount of serum).